



Health Professional Recommended  
Medicare Code: E0628

# Power Seat™

## Keeps People Independent at Home

**Uplift Power Seat™** is an electric-powered lift cushion that gently lifts and lowers users from an armchair or sofa.

### TARGET CUSTOMER

- Customers who need help getting up due to arthritis, Parkinson's, MS, ALS or other joint and neuromuscular conditions.
- Customers who do not want to replace their favorite armchair or sofa.
- Customers who own a lift chair but need something to use in other rooms.
- People who are mostly confined indoors - whether at home, visiting, or in assisted living.

### CUSTOMER BENEFITS

- Keeps people independent and at home longer.
- Can be stopped at any height needed.
- Conserves the users' energy for activities other than struggling to get up.

### BEST SEATS TO USE

- A favorite armchair
- The sofa - beside the arm
- Dining room chair
- Kitchen chair

### MEDICARE

- HCPCS Code: E0628, Medicare Part B
- Requires Certificate of Medical Necessity (CMN) from a physician.
- Conditions include arthritis, degenerative joint and neuromuscular disease.



*Bringing quality to life.*



# Power Seat™

## FEATURES

### EASY TO USE

The power lever acts like a switch. Simply click up to rise, push down to lower. The middle position will stop the lift at any time.



### PATENTED LevelLift™



The flexible seat curves to a flat surface that lifts users up on a level plane. This patented design provides for a safe lift up, with no forward dumping motion.

### POWER LEVERS



Power Seat comes with two detachable power levers that can be used on either the left or right side of the cushion.

### WATERPROOF



The breathable fabric cover is waterproof and machine washable in cold water. Comes in black.

### EASY CARE



The memory foam cushion has a water resistant coating that can be wiped down with mild soap and water.

### BUILT-IN HANDLE NON-SKID BOTTOM



The base has a built-in handle and a non-slip material covering that prevents the cushion from shifting on fabric seating surfaces.

### MEMORY FOAM CUSHION




Slow release memory foam provides additional comfort and prevention against pressure sores.



## Model PS1000

Lifts users up to 300lb

### UPLIFT POWER SEAT SPECIFICATIONS

<b>Model Number</b>	PS1000
<b>Medicare Code</b>	E0628
<b>Weight Capacity</b>	300lb/136kg
<b>Warranty</b>	1 Year Limited
<b>Weight</b>	Product 12lb/5.4kg Shipping 14lb/6.3kg
<b>Seat Dimensions</b>	Depth 19in/48cm Width 16in/41cm
<b>Seat Height</b>	Open 10in/25cm Flat 2in/5cm
<b>Electrical</b>	Power Supply 120 Vac 60Hz 0.6A CSA & C/US Approved Cord Length 10ft/3m 
<b>Patent Numbers</b>	5,316,370 & 6,702,383

**Uplift Technologies Inc.**  
#19-10 Morris Drive  
Dartmouth, NS Canada B3B 1K8  
[www.up-lift.com](http://www.up-lift.com)